

TUCSEN ISH500, ISH1000, TRUECHROME

Win 7/8/10 (64 bit)

i3 or better

4GB RAM or more

500GB HDD

USB 2.0, compatible with USB 3.0

GT Series

Win 7/8/10 (64 bit)

i3 or better

4GB RAM or more

500GB HDD

USB 2.0, compatible with USB 3.0

Michrome Series

Win 7/8/10 (64 bit)

i5 or better

8GB RAM or more

500GB HDD

USB 3.0

CARL ZEISS AXIOCAM

Cabinet : Minitower.Branded PC.

Intel® Core™ i5-4670 4-Core 3.0 GHz Intel®

Q87 chipset

4 GB DDR3-RAM

Graphics Adapter 1920 x 1200 Resolution, 32-bit true color, 1GB RAM, DirectX11.0 or higher -
OPTIONAL.HOWEVER GIVES BETTER RESULT IF HAVE THIS BOARD.

Monitor 20" TFT 1600 x 1200

Hard disk 1 TB SATA 3.0 Gbit/s

DVD writer

Microsoft® Windows® 7 SP1 x64 Professional or Ultimate or Microsoft® Windows® 10 x64 Enterprise
(not for Trigger board) No special customer adapted versions.

The user must be logged on as a member of the local "User" group.

For the software installation procedure local administrator rights are required.